



Umbrella	5pk	20pk	40pk
Retractor Medium	4870	4871	5162
Retractor Large	5256	5257	5258



**352 - OraSeal Syringe Kit**  
2 x 1.2 ml OraSeal Caulking syringes  
2 x 1.2 ml OraSeal Putty syringes  
4 x Black Mini tips  
20 x White Mac tips



1.2 ml syringe	4pk	20pk
Caulking	351	354
Putty	353	355



**311 - DermaDam Medium 0.20 mm 36pk**  
6"x 6" (15 cm x 15 cm)



**314 - DermaDam Heavy 0.25 mm 36pk**  
6"x 6" (15 cm x 15 cm)



**299 - DermaDam Medium Synthetic 0.20 mm 20pk**  
6"x 6" (15 cm x 15 cm)



**330 - DermaDam Medium Synthetic 0.20 mm 60pk**  
6"x 6" (15 cm x 15 cm)



**4100 - PropGard Kit**  
5 x Each PropGards large yellow and regular purple  
10 x Each Tongue guards left and right



PropGard	Regular 5pk	Large 5pk
Soft	4105	4103
Firm	4101	4104



**4102 - Tongue Guard 10pk**  
5 x Each right and left

Reliable Moisture Control,  
Clear Access, Superior Results



SCAN  
TO VIEW  
PRICING



## Umbrella™ Retractor

Developed to give you a clear treatment field while making patient comfort a top priority:

- 💧 Easy to place
- 💧 Disposable
- 💧 Naturally and gently helps the patient hold their mouth open without pulling or stretching their lips
- 💧 A new, innovative tongue-retraction design allows the tongue to comfortably rest behind the tongue guard, keeping it back and away from the working area
- 💧 Designed with anatomically placed/shaped bumpers, so clinicians can rest a hand on the patient's mouth without causing discomfort
- 💧 Provides relief and comfort to gaggers—it doesn't initiate the gag reflex for most
- 💧 Can be kept in place when checking bite
- 💧 Achieve a dry working environment



## PropGard™ Mouth Prop

Our PropGard mouth prop protects the tongue from trauma and prevents contamination from tongue and mouth closure.

- 💧 PropGard™ mouth prop protects adjacent teeth and tongue



## OraSeal™ Caulking and Putty

OraSeal™ Caulking and Putty materials effectively adhere to a dental dam, tissue, teeth, and metal even when under water or saliva.

- 💧 Excel at preventing leakage during treatment
- 💧 Protect the gingiva from exposure to peroxide or hydrofluoric acid
- 💧 Block out unwanted spaces for impressions
- 💧 In conjunction with the tear-resistant DermaDam rubber or DermaDam™ Synthetic dental dam, it's easy to achieve a dry working environment

## DermaDam™ Rubber Dam

DermaDam rubber dental dam is made from powder-free pure latex rubber. It has a low dermatitis potential and the lack of powder reduces the occurrence of allergic reactions.

- 💧 DermaDam Synthetic dental dam is not made with natural rubber latex but is designed to be just as flexible and durable

